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## Instruction Manual

### DDpolyX ML

### DDpolyX ML

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#### Material-related characteristics

##### Composition [wt.%]

Polymethylmethacrylat (PMMA)	> 99
Coloring pigments	< 1

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#### Explanation of symbols

Manufacturer	Date of manufacture
Consult instructions for use	Use by
<b>LOT</b>	<b>REF</b>
Batch code	Catalogue No.
Store dry	Protect from direct sun light
<b>R<sub>x</sub>only</b>	Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

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#### Instruction Manual

- Indications**  
DD polyX ML are tooth-coloured multilayer milling blanks made of PMMA for the manufacture of implant-supported temporary restorations, temporary crowns and bridges with up to two connected pontics in the front and side tooth area. Maximum recommended usage without check-up by dentist is 12 months.
- Contraindications**  
Parafunctional habits, intolerance to the contained components, model casting with clasps, bars, narrow sublingual bars and transversal connectors.
- Safety information**  
Please pay attention to the information in the safety data sheet in its current version. Avoid the inhalation of dust particles during processing. Wear protection gloves, safety glasses and a dust mask to prevent irritations of eyes, skin and respiratory system.
- Handling and storage**  
Only use clean, defect-free material out of undamaged packaging and make sure that the blanks are always stored in their original container in a dry and cool place, protected from light. Avoid contact with fluids.

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### 5. Designing

Handling this medical device should only be allowed to trained staff. When designing the objects, the minimum wall thickness should not be less than 0.8 mm in the cervical area and 1.5 mm in the occlusal area. The connection cross sections between crowns and pontics of a construction must not be less than 12 mm<sup>2</sup>. In general all constructions should be designed in a way to enable the compensation and the transfer of the mastication forces to the abutment teeth.

### 6. Milling and further processing

When positioning the blanks in the milling machine, pay attention to the correct colour gradient. The incisal layer is marked on the blank. For the design and positioning of the construction in the blank the individual layer heights can be taken from the following nesting table.

#### Nesting Table:

Only use single-edged milling tools in the milling machine or wet grinding procedures!

#### Further processing:

Tungsten carbide burs like the one you use for soft reline materials are suitable for rough adaption. For the finalization common cross-toothed milling tools can be used. Polishing should be carried out at low contact pressure, otherwise the resulting heat might lead to smearing or melting.

### 7. Cleaning

- ultrasonic bath (max. 40 °C) with mild cleaning agent

#### Do not use the following methods for cleaning:

- steam blasting
- autoclave method
- organic solvents and acids (e.g. acetone, hydrochloric acid)
- highly concentrated alcohols
- chemical denture cleaners
- abrasive cleaning agents

### 8. Mounting

#### Bonding extra-orally:

Choose primer suitable for the adhesive. Only use special adhesives, adapted to the material you want to combine with the PMMA. Pay attention to the instructions of the adhesive manufacturer.

#### Cementation:

For temporary cementations are recommended temporary cements or glass ionomer cements.

### 9. Disposal

Remains of PMMA blanks can be disposed of in the regular household garbage.

*We reserve the right to make changes because of the continuous development of our products. Please find the current version of the instruction manual on our website: [www.dentaldirekt.de](http://www.dentaldirekt.de).*

**This version replaces all previous versions.**