



VITA® A2

Liquid recommendation exemplified by the shade reference VITA® A2



DDBasicShade
Universal Dentine Liquid



DDProShadeZ
Special Dentine Liquid



DDProShadeC
Special Dentine Liquid



more grayish



more reddish



more grayish



more reddish



more grayish



more reddish



more grayish



more yellowish



Recommendation

Color concept: intensive chroma, higher gray value
Luminance: LESS brightness
DD recommendation: minimum layering (cut back), Full veneering, combination of both techniques

Color concept: less intensive chroma, warmer shade
Luminance: MORE brightness
DD recommendation: Further characterization: For monolithic crowns/bridges made of DD cubeX2® or DD cube ONE® an individual coloring (DD contrast® color and form concept) or minimal layering (cut back). For constructions made of DD Bio ZX2 or DD Bio Z a ceramic veneer.